## 顯微鏡生態探索 Microscope Eco Exploration

姓名 Name \_\_\_\_\_

組別 Group \_\_\_\_\_

日期 Date \_\_\_\_\_

### 學習目標 Learning goals

完成課程後，學生應能：After the course, students should be able to:

1. 有效收集水生微生物及動物樣本 Effectively collect aquatic microorganism and animal samples;
2. 製作裝載水生微生物的玻片 Prepare temporary mounts of aquatic microorganism samples;
3. 正確地使用複式顯微鏡及立體顯微鏡 Correctly use compound microscope and stereo microscope;
4. 明白光場、暗場和相位差不同之處 Understand the differences between bright field, dark field and phase contrast microscopy;
5. 辨別動物及認識其適應特徵 Identify animals found and understand their adaptation features;
6. 製作簡單科學報告 Make simple scientific report;
7. 欣賞大自然之美和尊重生物 Appreciate the wonder of nature and respect living things.

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In 2021, scientists from the Institute of Oceanography, Chinese Academy of Sciences, made a groundbreaking discovery by identifying marine microorganisms capable of effectively degrading plastic waste. They used electron microscopy to compare the degradation of plastic surfaces. The research findings have been successfully applied for a national invention patent, and it is hoped that they can contribute to solve the issue of marine plastic pollution in the future.

### 衣著 Clothing

1. 不應穿著短褲。穿著長袖上衣和長褲能更有效防止蚊蟲叮咬，亦可減低被太陽曬傷的機會。  
Shorts are not recommended. Long-sleeved shirt and trousers for better protection against mosquito and insect bites, as well as preventing sunburn.
2. 不應穿著拖鞋或涼鞋，而應穿著運動布鞋，以減低腳部受傷的機會。  
A pair of plimsolls for preventing injuries. Slippers and sandals are not recommended.

### 安全 Safety

1. 避免踏足陡斜和濕滑的岩石表面。  
Avoid stepping on steep and wet rock surfaces.
2. 切勿闖入植物生長茂密的地方，以免觸及蜂巢或被植物割傷。  
To prevent being attacked by wasps from disturbed wasp nests or being injured by plant leaves and thorns, do not get into places with dense vegetation.

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Search for the green or brown filamentous algae at the edge of the pond. Collect them by using a pair of forceps and a small vial.

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Use a dropper and a small vial to collect the sediment in the eco-pond.

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Use a dropper and a small vial to collect the sediment in the fish bucket.

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Collect two different species of tiny animals in Ho Koon trail or nearby area and place them in petri dishes.

## 實驗室工作：檢視樣本 Laboratory work: Sample examination

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Use a dropper to stir and then transfer 1 drop of water sample onto a slide. Then use a pair of forceps to collect a very small amount of algae and add into the water droplet. Cover with cover slip. Prepare at least two temporary mounts for observation.

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Use a dropper to suck up the sediment then expel most of the water until only sediments remain. Add 1 drop of water with sediments onto a slide. Cover with cover slip. Prepare at least two temporary mounts for observation.

#### ii. 複式顯微鏡檢視臨時載玻片 Observation of the temporary mounts slide under compound microscope

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Use a compound microscope to examine all prepared slides. Adjust the light intensity and the aperture diaphragm of the condenser to obtain the image with the best balance of brightness and contrast. Use bright field, dark field and phase contrast to observe the same sample. Identify the microorganisms observed. Note the shape, physical characteristics, behaviour, locomotion and other features of each species in the record sheet in page 3.

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Use a stereomicroscope to examine the collected animal samples. Adjust the light intensity to obtain the image with the best balance of brightness and contrast. Observe different parts of body (e.g. mouthpart, wings, eyes, leg, etc.) of the animal samples under different magnifications. Record all features observed in the record sheet in page 4.

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光場、暗場、相位差的影像比較 Comparison of the differences between bright field, dark field and phase contrast microscopy

	光場 Bright field	暗場 Dark field	相位差 Phase contrast
背景顏色 Background colour			
微生物結構在影像中的清晰度 Sharpness of the structures of the microorganisms in the image (0: 不能觀察 Cannot be observed 1: 能勉強觀察 Barely observed 2: 能清晰觀察 Clearly observed)	<input type="checkbox"/> 纖毛 Cilia <input type="checkbox"/> 鞭毛 Flagella <input type="checkbox"/> 體內組織 Body's tissue <input type="checkbox"/> 細胞器 Organelle	<input type="checkbox"/> 纖毛 Cilia <input type="checkbox"/> 鞭毛 Flagella <input type="checkbox"/> 體內組織 Body's tissue <input type="checkbox"/> 細胞器 Organelle	<input type="checkbox"/> 纖毛 Cilia <input type="checkbox"/> 鞭毛 Flagella <input type="checkbox"/> 體內組織 Body's tissue <input type="checkbox"/> 細胞器 Organelle

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## C. 分析和討論 Analysis and Discussion

- 簡略描述兩個水樣本所找到的微生物。  
Briefly describe the aquatic microorganisms found in the two water samples.
- 比較光場、暗場和相位差不同之處。  
Compare the differences between bright field, dark field and phase contrast microscopy.
- 繪畫一種在水樣本找到的原生動物高倍生物圖 ( X400 ) 。  
Draw a high power biological drawing (X400) of a protozoan that found in the water samples.
- 揀選一種你研究的動物樣本，介紹其外貌特徵與其行為習性的關係。  
Choose one of the animal sample investigated and introduce its external features in relation to its behaviour.